

Variable Angle EZ-SliceIR



This new accessory is very unique and convenient tool to cut the samples in several micron thickness. Angle can be set at any angle between 45 degree to 90 degree. Because of its unique mechanics, it can be used for not only soft samples like films but also hard samples like plastic.

Features:

- Sections of several micron thickness
- Adjustable cutting angle (45° degree - 90° degree)
- Cuts soft samples
- Cuts hard samples

ORDERING INFORMATION

Part Number

STJ-0205

Description

Variable Angle EZ-SliceIR